

GulfCare Water-in-Oil Reagent Determinate water content in lubricants/ marine liquid fuels

[Published in Mar 2023]



	Vater Reagent	re Water Reagent 20g x 4 tubes
Reagent A/Releaser/Cleaner 1 x 500mlPlastic Bottle 1 Unit10ml Syringe 2 UnitsGloves 4 Pairs50ml Graduated Tube 1 NumberMagnetic Pin, Steel Bal 1 Each	ttle Ige luated Tube	Bottle 1 Unityringe 2 Units 4 PairsGraduated Tube 1 Number

The GulfCare water reagent paste meets the necessary requirements for safe transportation and storage.

Each tube of paste weighs 20g and can be combined with 125ml of GulfCare Reagent A/Releaser/Cleaner to perform fourteen tests. With proper handling, one refill pack can facilitate up to 50 tests.

For optimal results, it is recommended to mix the paste with the liquid releaser immediately before conducting the test. To ensure consistency of results, it is advisable to avoid storing the prepared liquid for extended periods.

Additionally, the GulfCare water-in-oil test reagents can be utilized to determine water content in non-GulfCare brand test devices, provided that the same operational technique is applied.

Gulf Marine Pte. Ltd. e-mail : technical.engineer@gulf-marine.com web: www.gulf-marine.com